


## Recommended use

The Rifton Stander is a Class 1 medical device. It is intended to position people with disabilities in standing, providing the health and therapeutic benefits of supported standing. The Rifton Stander allows for growth and can be configured for prone or supine standing. The multi-position configuration allows both. A variety of padded supports, straps and adjustments allow comfortable and precise positioning for users with different physical needs. Configuration and initial setup should be done by a qualified medical professional.

## User and item dimensions

User dimensions – inches (cm)		Size 1
<b>Height</b>		<b>28 – 42 (71–107)*</b>
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate size. Choose the model that allows for growth. *The smallest users may require the S409 Sandal Lift Kit for best positioning.</p> <p><b>Important:</b> User's weight must not exceed the maximum working load.</p>		
Item dimensions – inches (cm)		Size 1
Base length		32½ (83)
Base width		21½ (55)
Horizontal body support height above floor		22–28 (56–71)
Sandal to knee height		6¾ – 11 (17 – 28)
Knee to hip height		6¼ – 10½ (16 – 27)
Sandal to top of headrest		31½ – 43 (80 – 109)
Distance between trunk laterals		5½ – 10 (14 – 25)
Distance between thigh medial and lateral support		2½ – 4 (6 – 10)
Basic item weight (Prone)		47 lbs (21 kg)
Basic item weight (Supine)		52 lbs (23.5 kg)
<b>Max. working load in lbs (kg)</b>		<b>75 (34 kg)</b>